

Easier to put on than elastic stockings,  
the CompreFit® is latex free, lightweight,  
durable and very comfortable.

**Options include:**

- BK (Below Knee)
- TC (Thigh Component)

*as seen below with BK (sold seperately)*

*BK includes CompreBoot™*

*TC includes a Hip Attachment  
and knee piece*



 **CompreBoot™**

 **CompreFit® -TC**  
UNIVERSAL

[www.BiaCare.com](http://www.BiaCare.com)

140 W. Washington Ave.  
Suite 100  
Zeeland, MI 49464

**BiaCare**  
Medical

**BiaCare**

Medical  
*Advanced Care for life*



 **CompreFit® -BK**  
UNIVERSAL

*Finally, compression wear that's*

*Easy*

*Effective*

*Affordable*

[www.BiaCare.com](http://www.BiaCare.com)

 **CompreFit®**

# CompreFit® UNIVERSAL

Most patients with Lymphedema, or other forms of edema, do not require custom compression garments. Designed specifically for those patients that have standard measurements, **BiaCare** Medical has designed standard compression wear to make life a little easier.

## Why CompreFit® may be the best choice for you...

- **Adjustable**  
*Means standard sizing*
- **Economical**  
*Saves time and money*
- **Effective**  
*Contains and reduces edema*
- **Easy Care**  
*Machine wash and dry*
- **Easy Use**  
*Means increased compliance*
- **Non-Elastic Compression**  
*Effective for day and nighttime*
- **Comfortable**  
*Made with soft fabric and neoprene*
- **Durable**  
*Lasts 2-3 times longer than bandages*



## CompreFit® -TC UNIVERSAL

The CompreFit® Universal, introduced in 2006, has experienced extremely positive reception from clinicians nationally. The unique properties of non-elastic compression allows the CompreFit® to provide a high working pressure (30-40 mmHg), while also providing a low resting pressure, making it ideal for day use and safe for night use. It is extremely durable, as well as easy to use and clean.

Many physicians prefer the CompreFit® to high compression elastic stockings for their patients.

## Adjustable

The innovative, Patent Pending design offers up to 4" (10cm) of circumferential adjustability. The length can also be adjusted up to 1" (2.5cm). Only a few measurements are needed to determine the correct size, with optional regular or tall lengths available.

The interlocking and overlapping strap design makes it easy to achieve consistent, gradient compression. The unique StrapKeeper™ feature aids garment handling and use.

## CompreFit® -BK UNIVERSAL

- Large, easy to use straps
- Straps overlap for "gauranteed gradient" compression
- Includes soft, 4-way stretch neoprene CompreBoot® and Cottonette stocking
- Offered in Black or Beige



## Economical

CompreFit® garments cost substantially less than other compression garments. considering the expected life of a CompreFit® is 2-3 times that of bandages, the cost of CompreFit® garments is lower than bandage wrapping. The CompreFit® saves time by eliminating wrapping/ re-rolling!

## CompreFit® -TC UNIVERSAL

- Includes soft, neoprene knee piece and Cottonette stocking
- Also includes Hip Attachment



BK Sold Separately

## Effective

CompreFit® garments are designed for daytime or nighttime use for individuals with the following indications:

- Lymphedema
- Post-surgical swelling
- Chronic Venous Insufficiency
- Phlebitis
- Other types of Edema

## Easy Care

CompreFit® garments are exceptionally durable and can be machine washed and dried (refer to care instructions).